

## 兔抗 MAPT(Ab-717/400) 多克隆抗体

中文名称：兔抗 MAPT(Ab-717/400) 多克隆抗体

英文名称： Anti-MAPT(Ab-717/400) rabbit polyclonal antibody

别名： TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; FTDP-17; PPP1R103

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： MAPT(Ab-717/400)

反应种属： Human Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes the microtubule-associated protein tau (MAPT) whose transcript undergoes complex, regulated alternative splicing, giving rise to several mRNA species. MAPT transcripts are differentially expressed in the nervous system, depending on stage of neuronal maturation and neuron type. MAPT gene mutations have been associated with several neurodegenerative disorders such as Alzheimer's disease, Pick's disease, frontotemporal dementia, cortico-basal degeneration and progressive supranuclear palsy.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MAPT(Ab-717/400)

<b>Immunogen:</b>	Synthesized non-phosphopeptide derived from human Tau around the phosphorylation site of serine 717/400 (V-V-S(p)-G-D).
<b>Full name:</b>	microtubule-associated protein tau
<b>Synonyms :</b>	TAU; MSTD; PPND; DDPAC; MAPTL; MTBT1; MTBT2; FTDP-17; PPP1R103
<b>SwissProt:</b>	P10636
<b>WB Predicted band size:</b>	79 kDa
<b>WB Positive control:</b>	Rat brain tissue lysate
<b>WB Recommended dilution:</b>	500-3000

