

## CAMK1G 抗原（重组蛋白）

中文名称：CAMK1G 抗原（重组蛋白）

英文名称：CAMK1G Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 23-277 amino acids of human CAMK1G

技术规格

<b>Full name:</b>	calcium/calmodulin-dependent protein kinase IG
<b>Synonyms:</b>	VWS1; CLICK3; CLICKIII; dJ272L16.1
<b>Swissprot:</b>	Q96NX5
<b>Gene Accession:</b>	BC032787
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known. Calcium/calmodulin-dependent protein kinase belonging to a proposed calcium-triggered signaling cascade. In vitro phosphorylates transcription factor CREB1 (By similarity).