

## 兔抗 MATN2 多克隆抗体

中文名称：兔抗 MATN2 多克隆抗体

英文名称：Anti-MATN2 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻 (-20℃)

宿 主：Rabbit

抗 原：MATN2

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the von Willebrand factor A domain containing protein family. This family of proteins is thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This protein contains five von Willebrand factor A domains. The specific function of this gene has not yet been determined. Two transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MATN2
<b>Immunogen:</b>	Synthetic peptide of human MATN2
<b>Full name:</b>	matrilin 2
<b>SwissProt:</b>	O00339

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human thyroid cancer and human colorectal cancer
<b>IHC Recommend dilution:</b>	25-100



