

## CASP2 抗原（重组蛋白）

中文名称：CASP2 抗原（重组蛋白）

英文名称：CASP2 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to N terminal 201 amino acids of human CASP2

技术规格

<b>Full name:</b>	Caspase-2
<b>Synonyms:</b>	ICH1, NEDD2, CASP-2, NEDD-2
<b>Swissprot:</b>	P42575
<b>Gene Accession:</b>	BC002427
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of the cysteine-aspartic acid protease (caspase) family. Caspases mediate cellular apoptosis through the proteolytic cleavage of specific protein substrates. The encoded protein may function in stress-induced cell death pathways, cell cycle maintenance, and the suppression of tumorigenesis. Increased expression of this gene may play a role in neurodegenerative disorders including Alzheimer's disease, Huntington's disease and temporal lobe epilepsy. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.