

CASP3 抗原（重组蛋白）

中文名称：CASP3 抗原（重组蛋白）

英文名称：CASP3 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to N terminal 269 amino acids of human CASP3

技术规格

Full name:	caspase 3, apoptosis-related cysteine peptidase
Synonyms:	CPP32, SCA-1, CPP32B
Swissprot:	P42574
Gene Accession:	BC016926
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein which is a member of the cysteine-aspartic acid protease (caspase) family. Sequential activation of caspases plays a central role in the execution-phase of cell apoptosis. Caspases exist as inactive proenzymes which undergo proteolytic processing at conserved aspartic residues to produce two subunits, large and small, that dimerize to form the active enzyme. This protein cleaves and activates caspases 6, 7 and 9, and the protein itself is processed by caspases 8, 9 and 10. It is the predominant caspase involved in the cleavage of amyloid-beta 4A precursor protein, which is associated with neuronal death in Alzheimer's disease. Alternative splicing

of this gene results in two transcript variants that encode the same protein.