

兔抗 MATN3 多克隆抗体

中文名称：兔抗 MATN3 多克隆抗体

英文名称： Anti-MATN3 rabbit polyclonal antibody

别名： HOA; OS2; EDM5; DIPOA; OADIP

抗原： MATN3

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of von Willebrand factor A domain containing protein family. This family of proteins is thought to be involved in the formation of filamentous networks in the extracellular matrices of various tissues. This protein contains two von Willebrand factor A domains; it is present in the cartilage extracellular matrix and has a role in the development and homeostasis of cartilage and bone. Mutations in this gene result in multiple epiphyseal dysplasia.
Applications:	ELISA, IHC
Name of antibody:	MATN3

Immunogen:	Fusion protein of human MATN3
Full name:	matrilin 3
Synonyms :	HOA; OS2; EDM5; DIPOA; OADIP
SwissProt:	O15232
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human cervical cancer and Human gastric cancer
IHC Recommend dilution:	25-100

