

## 兔抗 HHIPL2 多克隆抗体

中文名称：兔抗 HHIPL2 多克隆抗体

英文名称：Anti-HHIPL2 rabbit polyclonal antibody

别名：HHIP like 2; KIAA1822L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HHIPL2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The function of this protein remains unknown.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HHIPL2
<b>Immunogen:</b>	Fusion protein of human HHIPL2
<b>Full name:</b>	HHIP like 2
<b>Synonyms:</b>	KIAA1822L
<b>SwissProt:</b>	Q6UWX4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human lung cancer
<b>IHC Recommend dilution:</b>	50-200



