

## 兔抗 MAZ 多克隆抗体

中文名称：兔抗 MAZ 多克隆抗体

英文名称：Anti-MAZ rabbit polyclonal antibody

别名：PUR1; ZF87; Pur-1; SAF-1; SAF-2; SAF-3; Zif87; ZNF801

抗原：MAZ

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

### 技术规格

<b>Background:</b>	May function as a transcription factor with dual roles in transcription initiation and termination. Binds to two sites, ME1a1 and ME1a2, within the MYC promoter having greater affinity for the former. Also binds to multiple G/C-rich sites within the promoter of the Sp1 family of transcription factors. Regulates inflammation-induced expression of serum amyloid A proteins.
<b>Applications:</b>	WB
<b>Name of antibody:</b>	MAZ
<b>Immunogen:</b>	Synthesized peptide derived from internal of human MAZ.
<b>Full name:</b>	MYC-associated zinc finger protein (purine-binding transcr

	iption factor)
<b>Synonyms :</b>	PUR1; ZF87; Pur-1; SAF-1; SAF-2; SAF-3; Zif87; ZNF801
<b>SwissProt:</b>	P56270
<b>WB Predicted band size:</b>	49 kDa
<b>WB Positive control:</b>	K562 cells lysate
<b>WB Recommended dilution:</b>	500-3000

