

兔抗 MBL2 多克隆抗体

中文名称：兔抗 MBL2 多克隆抗体

英文名称：Anti-MBL2 rabbit polyclonal antibody

别名：MBL; MBP; MBP1; MBPD; MBL2D; MBP-C; COLEC1; HSMBPC

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：MBL2

反应种属：Human, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes the soluble mannose-binding lectin or mannose-binding protein found in serum. The protein encoded belongs to the collectin family and is an important element in the innate immune system. The protein recognizes mannose and N-acetylglucosamine on many microorganisms, and is capable of activating the classical complement pathway. Deficiencies of this gene have been associated with susceptibility to autoimmune and infectious diseases.
Applications:	ELISA, IHC
Name of antibody:	MBL2

Immunogen:	Synthetic peptide of human MBL2
Full name:	mannose binding lectin 2
Synonyms:	MBL; MBP; MBP1; MBPD; MBL2D; MBP-C; COLEC1; HSMBP C
SwissProt:	P11226
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and Human esophagus cancer
IHC Recommend dilution:	50-100



