

兔抗 MBNL3 多克隆抗体

中文名称：兔抗 MBNL3 多克隆抗体

英文名称： Anti-MBNL3 rabbit polyclonal antibody

别名： CHCR; MBLX; MBXL; MBLX39

相关类别： 一抗

抗原： MBNL3

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the muscleblind-like family of proteins. The encoded protein may function in regulation of alternative splicing and may play a role in the pathophysiology of myotonic dystrophy. Alternatively spliced transcript variants have been described.
Applications:	ELISA, IHC
Name of antibody:	MBNL3
Immunogen:	Fusion protein of human MBNL3
Full name:	muscleblind-like splicing regulator 3
Synonyms :	CHCR; MBLX; MBXL; MBLX39
SwissProt:	Q9NUK0

ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and human cervical cancer
IHC Recommend dilution:	50-200

