

## 兔抗 MBOAT4 多克隆抗体

中文名称: 兔抗 MBOAT4 多克隆抗体

英文名称: Anti-MBOAT4 rabbit polyclonal antibody

别名: membrane bound O-acyltransferase domain containing 4; GOAT; OACT4; FKSG89

抗原: MBOAT4

储存: 冷冻 (-20°C)

抗原: MBOAT4

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Mediates the octanoylation of ghrelin at 'Ser-3'. Can use a variety of fatty acids as substrates including octanoic acid, decanoic acid and tetradecanoic acid.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MBOAT4
<b>Immunogen:</b>	Synthetic peptide of human MBOAT4
<b>Full name:</b>	membrane bound O-acyltransferase domain containing 4
<b>Synonyms:</b>	GOAT; OACT4; FKSG89
<b>SwissProt:</b>	Q96T53

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human colorectal cancer
<b>IHC Recommend dilution:</b>	30-150



