

兔抗 HIPK1 多克隆抗体

中文名称：兔抗 HIPK1 多克隆抗体

英文名称：Anti-HIPK1 rabbit polyclonal antibody

别名：Myak; Nbak2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HIPK1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the Ser/Thr family of protein kinases and HIPK subfamily. It phosphorylates homeodomain transcription factors and may also function as a co-repressor for homeodomain transcription factors. Alternative splicing results in four transcript variants encoding four distinct isoforms.
Applications:	ELISA, IHC
Name of antibody:	HIPK1
Immunogen:	Synthetic peptide of human HIPK1
Full name:	homeodomain interacting protein kinase 1
Synonyms :	Myak; Nbak2

SwissProt:	Q86Z02
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	50-100

