

## 兔抗 MBOAT7 多克隆抗体

中文名称：兔抗 MBOAT7 多克隆抗体

英文名称： Anti-MBOAT7 rabbit polyclonal antibody

别名： membrane bound O-acyltransferase domain containing 7; BB1; LRC4; LENG4; LPIAT; MBOA7; MRT57; OACT7; hMBOA-7

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MBOAT7

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the membrane-bound O-acyltransferases family of integral membrane proteins that have acyltransferase activity. The encoded protein is a lysophosphatidylinositol acyltransferase that has specificity for arachidonoyl-CoA as an acyl donor. This protein is involved in the reacylation of phospholipids as part of the phospholipid remodeling pathway known as the Land cycle. Alternative splicing results in multiple transcript variants.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	MBOAT7
<b>Immunogen:</b>	Fusion protein of human MBOAT7
<b>Full name:</b>	membrane bound O-acyltransferase domain containing 7
<b>Synonyms:</b>	BB1; LRC4; LENG4; LPIAT; MBOA7; MRT57; OACT7; hMBOA-7
<b>SwissProt:</b>	Q96N66
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer and Human tonsil
<b>IHC Recommend dilution:</b>	10-50



