

兔抗 MBTPS1 多克隆抗体

中文名称：兔抗 MBTPS1 多克隆抗体

英文名称： Anti-MBTPS1 rabbit polyclonal antibody

别名： S1P; PCSK8; SKI-1

抗原： MBTPS1

储存： 冷冻（-20℃） 避光

宿主： Rabbit

反应种属： Human Mouse Rat

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene encodes a member of the subtilisin-like propeptidase family, which includes proteases that process protein and peptide precursors trafficking through regulated or constitutive branches of the secretory pathway. The encoded protein undergoes an initial autocatalytic processing event in the ER to generate a heterodimer which exits the ER and sorts to the cis/medial-Golgi where a second autocatalytic event takes place and the catalytic activity is acquired.
Applications:	WB
Name of antibody:	MBTPS1

Immunogen:	Synthesized peptide derived from C-terminal of human MB TPS1.
Full name:	membrane-bound transcription factor peptidase, site 1
Synonyms :	S1P; PCSK8; SKI-1
SwissProt:	Q14703
WB Predicted band size:	118 kDa
WB Positive control:	Jurkat and COLO205 lysates
WB Recommended dilution:	500-3000

