

## 兔抗 HIPK4 多克隆抗体

中文名称：兔抗 HIPK4 多克隆抗体

英文名称：Anti-HIPK4 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

抗 原：HIPK4

宿 主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	HIPK4 (Homeodomain-interacting protein kinase 4) is a 616 amino acid protein that localizes to the cytoplasm and exists as two isoforms that are produced by alternative splicing events. Like other members of the HIPK family, HIPK4 functions as a protein kinase that catalyzes the ATP-dependent phosphorylation of target proteins and is thought to act as a corepressor of target transcription factors.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HIPK4
<b>Immunogen:</b>	Fusion protein of human HIPK4
<b>Full name:</b>	homeodomain interacting protein kinase 4
<b>SwissProt:</b>	Q8NE63
<b>ELISA Recommended dilutio</b>	2000-5000

<b>n:</b>	
<b>IHC positive control:</b>	Human thyroid cancer
<b>IHC Recommend dilution:</b>	25-100

