

兔抗 FUT4 多克隆抗体

中文名称：兔抗 FUT4 多克隆抗体

英文名称： Anti-FUT4 rabbit polyclonal antibody

别名： fucosyltransferase 4; LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV

相关类别： 一抗

抗原： FUT4

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The product of this gene transfers fucose to N-acetyllactosamine polysaccharides to generate fucosylated carbohydrate structures. It catalyzes the synthesis of the non-sialylated antigen, Lewis x (CD15).
Applications:	ELISA, IHC
Name of antibody:	FUT4
Immunogen:	Synthetic peptide of human FUT4
Full name:	fucosyltransferase 4
Synonyms:	LeX; CD15; ELFT; FCT3A; FUTIV; SSEA-1; FUC-TIV
SwissProt:	P22083
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human colorectal cancer

IHC Recommend dilution: 30-150





