

兔抗 MC4R 多克隆抗体

中文名称：兔抗 MC4R 多克隆抗体

英文名称： Anti-MC4R rabbit polyclonal antibody

相关类别： 一抗

储 存： 冷冻 (-20℃)

宿 主： Rabbit

抗 原： MC4R

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a membrane-bound receptor and member of the melanocortin receptor family. The encoded protein interacts with adrenocorticotrophic and MSH hormones and is mediated by G proteins. This is an intronless gene. Defects in this gene are a cause of autosomal dominant obesity
Applications:	ELISA, IHC
Name of antibody:	MC4R
Immunogen:	Synthetic peptide of human MC4R
Full name:	Melanocortin 4 receptor
SwissProt:	P32245
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human ovarian cancer and Human lung cancer

IHC Recommend dilution: 50-200

