

兔抗 FUT10 多克隆抗体

中文名称：兔抗 FUT10 多克隆抗体

英文名称：Anti-FUT10 rabbit polyclonal antibody

别名：fucosyltransferase 10; FUCTX

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FUT10

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|------------------------------------|----------------------------------|
| Background: | Probable fucosyltransferase. |
| Applications: | ELISA, WB, IHC |
| Name of antibody: | FUT10 |
| Immunogen: | Synthetic peptide of human FUT10 |
| Full name: | fucosyltransferase 10 |
| Synonyms: | FUCTX |
| SwissProt: | Q6P4F1 |
| ELISA Recommended dilution: | 5000-10000 |
| IHC positive control: | Human esophagus cancer |
| IHC Recommend dilution: | 40-200 |
| WB Predicted band size: | 56 kDa |
| WB Positive control: | Mouse liver tissue lysate |

WB Recommended dilution:

500-2000



