

兔抗 MCCC1 多克隆抗体

中文名称: 兔抗 MCCC1 多克隆抗体

英文名称: Anti-MCCC1 rabbit polyclonal antibody

别名: methylcrotonoyl-CoA carboxylase 1; MCCA; MCC-B

抗原: MCCC1

储存: 冷冻 (-20°C)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes the large subunit of 3-methylcrotonyl-CoA carboxylase. This enzyme functions as a heterodimer and catalyzes the carboxylation of 3-methylcrotonyl-CoA to form 3-methylglutaconyl-CoA. Mutations in this gene are associated with 3-Methylcrotonylglycinuria, an autosomal recessive disorder of leucine catabolism.
Applications:	ELISA, WB, IHC
Name of antibody:	MCCC1
Immunogen:	Fusion protein of human MCCC1
Full name:	methylcrotonoyl-CoA carboxylase 1
Synonyms:	MCCA; MCC-B

SwissProt:	Q96RQ3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	40-200
WB Predicted band size:	80 kDa
WB Positive control:	293T, HepG2, A172 cell, Mouse heart tissue, Jurkat cell lysates
WB Recommended dilution:	500-2000



