

## 兔抗 MCHR2 多克隆抗体

中文名称: 兔抗 MCHR2 多克隆抗体

英文名称: Anti-MCHR2 rabbit polyclonal antibody

别名: melanin concentrating hormone receptor 2; SLT; MCH2; MCH2R; GPR145; GPRv17; MCH-2R; MCH-R2; MCHR-2

抗原: MCHR2

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

|                                    |  |
|------------------------------------|--|
| <b>Background:</b>                 | Receptor for melanin-concentrating hormone, coupled to G proteins that activate phosphoinositide hydrolysis. |
| <b>Applications:</b>               | ELISA, IHC   |
| <b>Name of antibody:</b>           | MCHR2  |
| <b>Immunogen:</b>                  | Synthetic peptide of human MCHR2   |
| <b>Full name:</b>                  | melanin concentrating hormone receptor 2   |
| <b>Synonyms:</b>                   | SLT; MCH2; MCH2R; GPR145; GPRv17; MCH-2R; MCH-R2; MCHR-2   |
| <b>SwissProt:</b>                  | Q969V1   |
| <b>ELISA Recommended dilution:</b> | 5000-10000   |

|                                |                      |
|--------------------------------|----------------------|
| <b>IHC positive control:</b>   | Human thyroid cancer |
| <b>IHC Recommend dilution:</b> | 20-100               |

