

CEACAM21 抗原（重组蛋白）

中文名称：CEACAM21 抗原（重组蛋白）

英文名称：CEACAM21 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 35-191 amino acids of human CEACAM21

技术规格

Full name:	carcinoembryonic antigen-related cell adhesion molecule 21
Synonyms:	CEACAM3, R29124_1
Swissprot:	Q3KPI0
Gene Accession:	BC012001
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	CEACAM21 is a gene which encodes carcinoembryonic antigen. It belongs to the immunoglobulin superfamily-CEA family. Three isoforms of the human protein are produced by alternative splicing. High levels of expression in hematopoietic system, brain and muscle.