

## CFLAR 抗原（重组蛋白）

中文名称：CFLAR 抗原（重组蛋白）

英文名称：CFLAR Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to N terminal 200 amino acids of human CFLAR

技术规格

<b>Full name:</b>	CASP8 and FADD-like apoptosis regulator
<b>Synonyms:</b>	CASH; FLIP; MRIT; CLARP; FLAME; Casper; FLAME1; c-FLIP; FLAME-1 ; I-FLICE; c-FLIPL; c-FLIPR; c-FLIPS; CASP8AP1
<b>Swissprot:</b>	O15519
<b>Gene Accession:</b>	BC001602
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a regulator of apoptosis and is structurally similar to caspase-8. However, the encoded protein lacks caspase activity and appears to be itself cleaved into two peptides by caspase-8. Several transcript variants encoding different isoforms have been found for this gene, and partial evidence for several more variants exists.