

兔抗 FXR1 多克隆抗体

中文名称：兔抗 FXR1 多克隆抗体

英文名称：Anti-FXR1 rabbit polyclonal antibody

别名：FXR1P

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：FXR1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

| | |
|--------------------------|---|
| Background: | The protein encoded by this gene is an RNA binding protein that interacts with the functionally-similar proteins FMR1 and FXR2. These proteins shuttle between the nucleus and cytoplasm and associate with polyribosomes, predominantly with the 60 S ribosomal subunit. Three transcript variants encoding different isoforms have been found for this gene |
| Applications: | ELISA, IHC |
| Name of antibody: | FXR1 |
| Immunogen: | Synthetic peptide of human FXR1 |

| | |
|------------------------------------|---|
| Full name: | fragile X mental retardation, autosomal homolog 1 |
| Synonyms : | FXR1P |
| SwissProt: | P51114 |
| ELISA Recommended dilution: | 1000-2000 |
| IHC positive control: | Human thyroid cancer and Human gastric cancer |
| IHC Recommend dilution: | 25-100 |

