

## 兔抗 FXYD6 多克隆抗体

中文名称：兔抗 FXYD6 多克隆抗体

英文名称：Anti-FXYD6 rabbit polyclonal antibody

别名：FXYD domain containing ion transport regulator 6

相关类别：一抗

储存：冷冻（-20℃）

抗原：FXYD6

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the FXYD family of transmembrane proteins. This particular protein encodes phosphohippolin, which likely affects the activity of Na,K-ATPase. Multiple alternatively spliced transcript variants encoding the same protein have been described. Related pseudogenes have been identified on chromosomes 10 and X. Read-through transcripts have been observed between this locus and the downstream sodium/potassium-transporting ATPase subunit gamma (FXYD2, GeneID 486) locus.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	FXYD6
<b>Immunogen:</b>	Fusion protein of human FXYD6
<b>Full name:</b>	FXYD domain containing ion transport regulator 6

<b>SwissProt:</b>	Q9H0Q3
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	11 kDa
<b>WB Positive control:</b>	Human cerebrum tissue, Mouse brain tissue, Human cereb ella tissue lysates
<b>WB Recommended dilution:</b>	500-2000



