

兔抗 HK3 多克隆抗体

- 中文名称： 兔抗 HK3 多克隆抗体
- 英文名称： Anti-HK3 rabbit polyclonal antibody
- 别名： HXK3; HKIII
- 相关类别： 一抗
- 抗原： HK3
- 储存： 冷冻（-20℃）
- 宿主： Rabbit
- 反应种属： Human, Mouse, Rat
- 标记物： Unconjugate
- 克隆类型： rabbit polyclonal

技术规格

Background:	Hexokinases phosphorylate glucose to produce glucose-6-phosphate, the first step in most glucose metabolism pathways. This gene encodes hexokinase 3. Similar to hexokinases 1 and 2, this allosteric enzyme is inhibited by its product glucose-6-phosphate.
Applications:	ELISA, IHC
Name of antibody:	HK3
Immunogen:	Fusion protein of human HK3
Full name:	hexokinase 3 (white cell)
Synonyms :	HXK3; HKIII
SwissProt:	P52790

ELISA Recommended dilution:	1000-5000
IHC positive control:	Human cervical cancer and human ovarian cancer
IHC Recommend dilution:	50-200

