

## CLEC1B 抗原（重组蛋白）

中文名称：CLEC1B 抗原（重组蛋白）

英文名称：CLEC1B Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-196 amino acids of human CLEC1B

技术规格

<b>Full name:</b>	C-type lectin domain family 1 member B
<b>Synonyms:</b>	CLEC2; CLEC2B; PRO1384; QDED721; 1810061I13Rik
<b>Swissprot:</b>	Q9P126
<b>Gene Accession:</b>	BC029554
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Natural killer (NK) cells express multiple calcium-dependent (C-type) lectin-like receptors, such as CD94 (KLRD1; MIM 602894) and NKG2D (KLRC4; MIM 602893), that interact with major histocompatibility complex class I molecules and either inhibit or activate cytotoxicity and cytokine secretion. CLEC2 is a C-type lectin-like receptor expressed in myeloid cells and NK cells (Colonna et al., 2000 [PubMed 10671229]).