

兔抗 FYCO1 多克隆抗体

中文名称：兔抗 FYCO1 多克隆抗体

英文名称： Anti-FYCO1 rabbit polyclonal antibody

别名： FYVE and coiled-coil domain containing 1; CATC2; RUFY3; ZFYVE7; CTRCT18

抗原： FYCO1

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a protein that contains a RUN domain, FYVE-type zinc finger domain and Golgi dynamics (GOLD) domain. The encoded protein plays a role in microtubule plus end-directed transport of autophagic vesicles through interactions with the small GTPase Rab7, phosphatidylinositol-3-phosphate (PI3P) and the autophagosome marker LC3. Mutations in this gene are a cause of autosomal recessive congenital cataract-2 (CATC2).
Applications:	ELISA, IHC
Name of antibody:	FYCO1

Immunogen:	Fusion protein of human FYCO1
Full name:	FYVE and coiled-coil domain containing 1
Synonyms:	CATC2; RUFY3; ZFYVE7; CTRCT18
SwissProt:	Q9BQS8
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human cervical cancer
IHC Recommend dilution:	100-200



