

CLEC4A 抗原（重组蛋白）

中文名称：CLEC4A 抗原（重组蛋白）

英文名称：CLEC4A Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 70-237 amino acids of human CLEC4A

技术规格

Full name:	C-type lectin domain family 4, member A
Synonyms:	DCIR, LLIR, DDB27, CLECSF6, HDCGC13P
Swissprot:	Q9UMR7
Gene Accession:	BC117439
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. The encoded type 2 transmembrane protein may play a role in inflammatory and immune response. Multiple transcript variants encoding distinct isoforms have been identified for this gene. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region.

mlbio 梅联生物
Good elisakit producers