

兔抗 MCM5 多克隆抗体

中文名称：兔抗 MCM5 多克隆抗体

英文名称： Anti-MCM5 rabbit polyclonal antibody

别名： CDC46; P1-CDC46

抗原： MCM5

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is structurally very similar to the CDC46 protein from <i>S. cerevisiae</i> , a protein involved in the initiation of DNA replication. The encoded protein is a member of the MCM family of chromatin-binding proteins and can interact with at least two other members of this family. The encoded protein is upregulated in the transition from the G0 to G1/S phase of the cell cycle and may actively participate in cell cycle regulation.
Applications:	ELISA, WB, IHC
Name of antibody:	MCM5
Immunogen:	Fusion protein of human MCM5

Full name:	minichromosome maintenance complex component 5
Synonyms:	CDC46; P1-CDC46
SwissProt:	P33992
IHC positive control:	Human gastric cancer and human liver cancer
IHC Recommend dilution:	50-200
WB Predicted band size:	82 kDa
WB Positive control:	HepG2, K562, Jurkat, 231 and hela cells
WB Recommended dilution:	500-2000



