

CLEC4D 抗原（重组蛋白）

中文名称：CLEC4D 抗原（重组蛋白）

英文名称：CLEC4D Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 39-215 amino acids of human CLEC4D

技术规格

Full name:	C-type lectin domain family 4 member D
Synonyms:	MCL; MPCL; CD368; CLEC6; CLEC-6; CLECSF8; Dectin-3
Swissprot:	Q8WXI8
Gene Accession:	BC032313
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. Members of this family share a common protein fold and have diverse functions, such as cell adhesion, cell-cell signalling, glycoprotein turnover, and roles in inflammation and immune response. This gene is closely linked to other CTL/CTLD superfamily members on chromosome 12p13 in the natural killer gene complex region.