

CLEC7A 抗原（重组蛋白）

中文名称：CLEC7A 抗原（重组蛋白）

英文名称：CLEC7A Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1-77 amino acids of human CLEC7A

技术规格

Full name:	C-type lectin domain containing 7A
Synonyms:	BGR; CD369; CANDF4; SCARE2; DECTIN1; CLECSF12
Swissprot:	Q9BXN2
Gene Accession:	BC013385
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a member of the C-type lectin/C-type lectin-like domain (CTL/CTLD) superfamily. The encoded glycoprotein is a small type II membrane receptor with an extracellular C-type lectin-like domain fold and a cytoplasmic domain with an immunoreceptor tyrosine-based activation motif. It functions as a pattern-recognition receptor that recognizes a variety of beta-1,3-linked and beta-1,6-linked glucans from fungi and plants, and in this way plays a role in innate immune response. Alternate transcriptional splice variants, encoding different isoforms, have been characterized. This gene is closely linked to other CTL/CTLD superfamily members.

rs on chromosome 12p13 in the natural killer gene complex region.