

## CMTM2 抗原（重组蛋白）

中文名称：CMTM2 抗原（重组蛋白）

英文名称：CMTM2 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to N terminal 115 amino acids of human CMTM2

技术规格

<b>Full name:</b>	CKLF-like MARVEL transmembrane domain containing 2
<b>Synonyms:</b>	CKLFSF2
<b>Swissprot:</b>	Q8TAZ6
<b>Gene Accession:</b>	BC025354
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene belongs to the chemokine-like factor gene superfamily, a novel family that links the chemokine and the transmembrane 4 su perfamilies of signaling molecules. The protein encoded by this gen e may play an important role in testicular development.