

## 兔抗 FYTTD1 多克隆抗体

中文名称: 兔抗 FYTTD1 多克隆抗体

英文名称: Anti-FYTTD1 rabbit polyclonal antibody

别名: forty-two-three domain containing 1; UIF

相关类别: 一抗

储存: 冷冻 (-20℃)

宿主: Rabbit

抗原: FYTTD1

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	Required for mRNA export from the nucleus to the cytoplasm. Acts as an adapter that uses the DDX39B/UAP56-NFX1 pathway to ensure efficient mRNA export and delivering to the nuclear pore. Associates with spliced and unspliced mRNAs simultaneously with ALYREF/THOC4.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	FYTTD1
<b>Immunogen:</b>	Fusion protein of human FYTTD1
<b>Full name:</b>	forty-two-three domain containing 1
<b>Synonyms:</b>	UIF
<b>SwissProt:</b>	Q96QD9

<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	30-150



