

## COL2A1 抗原（重组蛋白）

中文名称：COL2A1 抗原（重组蛋白）

英文名称：COL2A1 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 1220-1487 amino acids of human COL2A1

技术规格

<b>Full name:</b>	collagen, type II, alpha 1
<b>Synonyms:</b>	AOM, ANFH, SEDC, COL11A3
<b>Swissprot:</b>	P02458
<b>Gene Accession:</b>	BC007252
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes the alpha-1 chain of type II collagen, a fibrillar collagen found in cartilage and the vitreous humor of the eye. Mutations in this gene are associated with achondrogenesis, chondrodysplasia, early onset familial osteoarthritis, SED congenita, Langer-Saldino achondrogenesis, Kniest dysplasia, Stickler syndrome type I, and spondyloepimetaphyseal dysplasia Strudwick type. In addition, defects in processing chondrocalcin, a calcium binding protein that is the C-propeptide of this collagen molecule, are also associated with chondrodysplasia. There are two transcripts identified for this gene.

**mlbio** 梅联生物  
Good elisakit producers