

## CRMP1 抗原（重组蛋白）

中文名称：CRMP1 抗原（重组蛋白）

英文名称：CRMP1 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to C terminal 249 amino acids of human CRMP1

技术规格

<b>Full name:</b>	collapsin response mediator protein 1
<b>Synonyms:</b>	DRP1, DRP-1, CRMP-1, DPYSL1, ULIP-3
<b>Swissprot:</b>	Q14194
<b>Gene Accession:</b>	BC000252
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	This gene encodes a member of a family of cytosolic phosphoproteins expressed exclusively in the nervous system. The encoded protein is thought to be a part of the semaphorin signal transduction pathway implicated in semaphorin-induced growth cone collapse during neural development. Alternative splicing results in multiple transcript variants.