

兔抗 FZD2 多克隆抗体

中文名称：兔抗 FZD2 多克隆抗体

英文名称： Anti-FZD2 rabbit polyclonal antibody

别名： Fz2; fz-2; fzE2; hFz2

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： FZD2

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This intronless gene is a member of the frizzled gene family. Members of this family encode seven-transmembrane domain proteins that are receptors for the wingless type MMTV integration site family of signaling proteins. This gene encodes a protein that is coupled to the beta-catenin canonical signaling pathway. Competition between the wingless-type MMTV integration site family, member 3A and wingless-type MMTV integration site family, member 5A gene products for binding of this protein is thought to regulate the beta-catenin-dependent and -independent pathways.

Applications:	ELISA, IHC
Name of antibody:	FZD2
Immunogen:	Synthetic peptide of human FZD2
Full name:	led family receptor 2
Synonyms :	Fz2; fz-2; fzE2; hFz2
SwissProt:	Q14332
ELISA Recommended dilution:	1000-5000
IHC positive control:	Human renal cancer
IHC Recommend dilution:	25-100

