

兔抗 MED13 多克隆抗体

中文名称: 兔抗 MED13 多克隆抗体

英文名称: Anti-MED13 rabbit polyclonal antibody

别 名: ARC250; THRAP1; DRIP250; HSPC221; TRAP240

抗 原: MED13

储 存: 冷冻(-20℃)

宿 主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a component of the mediator comple x (also known as TRAP, SMCC, DRIP, or ARC), a transcrip tional coactivator complex thought to be required for the expression of almost all genes. The mediator complex is recruited by transcriptional activators or nuclear recept ors to induce gene expression, possibly by interacting with RNA polymerase II and promoting the formation of a transcriptional pre-initiation complex.
Applications:	ELISA, WB
Name of antibody:	MED13
Immunogen:	Synthetic peptide of human MED13



Full name:	mediator complex subunit 13
Synonyms:	ARC250; THRAP1; DRIP250; HSPC221; TRAP240
SwissProt:	Q9UHV7
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	239 kDa
WB Positive control:	K562 cells
WB Recommended dilution:	500-2000

