

DKK1 抗原（重组蛋白）

中文名称： DKK1 抗原（重组蛋白）

英文名称： DKK1 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 32-266 amino acids of human DKK1

技术规格

Full name:	dickkopf WNT signaling pathway inhibitor 1
Synonyms:	SK; DKK-1
Swissprot:	O94907
Gene Accession:	BC001539
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that is a member of the dickkopf family. It is a secreted protein with two cysteine rich regions and is involved in embryonic development through its inhibition of the WNT signaling pathway. Elevated levels of DKK1 in bone marrow plasma and peripheral blood is associated with the presence of osteolytic bone lesions in patients with multiple myeloma.