

兔抗 G3BP2 多克隆抗体

中文名称：兔抗 G3BP2 多克隆抗体

英文名称： Anti-G3BP2 rabbit polyclonal antibody

储 存： 冷冻（-20℃） 避光

抗 原： G3BP2

宿 主： Rabbit

反应种属： Human Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	Probable scaffold protein that may be involved in mRNA transport Potential.
Applications:	WB
Name of antibody:	G3BP2
Immunogen:	Synthesized peptide derived from C-terminal of human G3BP2.
Full name:	GTPase activating protein (SH3 domain) binding protein 2
Synonyms :	0
SwissProt:	Q9UN86
WB Predicted band size:	54 kDa
WB Positive control:	A549 cells, Huvec cells and Hela cells lysates
WB Recommended dilution:	500-3000

