

## 兔抗 HLTF 多克隆抗体

中文名称：兔抗 HLTF 多克隆抗体

英文名称：Anti-HLTF rabbit polyclonal antibody

别名：ZBU1; HLTF1; RNF80; HIP116; SNF2L3; HIP116A; SMARCA3

相关类别：一抗

储存：冷冻（-20℃）

抗原：HLTF

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the SWI/SNF family. Members of this family have helicase and ATPase activities and are thought to regulate transcription of certain genes by altering the chromatin structure around those genes. The encoded protein contains a RING finger DNA binding motif. Two transcript variants encoding the same protein have been found for this gene. However, use of an alternative translation start site produces an isoform that is truncated at the N-terminus compared to the full-length protein.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	HLTF

<b>Immunogen:</b>	Synthetic peptide of human HLTF
<b>Full name:</b>	helicase-like transcription factor
<b>Synonyms :</b>	ZBU1; HLTF1; RNF80; HIP116; SNF2L3; HIP116A; SMARCA3
<b>SwissProt:</b>	Q14527
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human thyroid cancer and Human gastric cancer
<b>IHC Recommend dilution:</b>	50-100

