

DKK4 抗原（重组蛋白）

中文名称： DKK4 抗原（重组蛋白）

英文名称： DKK4 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 21-224 amino acids of human DKK4

技术规格

Full name:	dickkopf WNT signaling pathway inhibitor 4
Synonyms:	DKK-4
Swissprot:	Q9UBT3
Gene Accession:	BC107047
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	This gene encodes a protein that is a member of the dickkopf family. The secreted protein contains two cysteine rich regions and is involved in embryonic development through its interactions with the Wnt signaling pathway. Activity of this protein is modulated by binding to the Wnt co-receptor and the co-factor kremen 2.