

DNMT3L 抗原（重组蛋白）

中文名称：DNMT3L 抗原（重组蛋白）

英文名称：DNMT3L Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to a region derived from 87-386 amino acids of human DNMT3L

技术规格

Full name:	DNA (cytosine-5-)-methyltransferase 3-like
Swissprot:	Q9UJW3
Gene Accession:	BC002560
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	CpG methylation is an epigenetic modification that is important for embryonic development, imprinting, and X-chromosome inactivation . Studies in mice have demonstrated that DNA methylation is required for mammalian development. This gene encodes a nuclear protein with similarity to DNA methyltransferases, but is not thought to function as a DNA methyltransferase as it does not contain the amino acid residues necessary for methyltransferase activity. However, it does stimulate de novo methylation by DNA cytosine methyltransferase 3 alpha and is thought to be required for the establishment of maternal genomic imprints. This protein also mediates transcriptional repression through interaction with histone deacetylase

1. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.