

兔抗 MED18 多克隆抗体

中文名称：兔抗 MED18 多克隆抗体

英文名称： Anti-MED18 rabbit polyclonal antibody

别名： p28b

储存： 冷冻（-20℃）

抗原： MED18

宿主： Rabbit

反应种属： Human, Mouse

相关类别： 一抗

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	MED18 is a component of the Mediator complex, which is a coactivator for DNA-binding factors that activate transcription via RNA polymerase II (Sato et al., 2003 [PubMed 12584197]).
Applications:	ELISA, IHC
Name of antibody:	MED18
Immunogen:	Full length fusion protein
Full name:	mediator complex subunit 18
Synonyms :	p28b
SwissProt:	Q9BUE0
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human liver cancer and human thyroid cancer
IHC Recommend dilution:	20-100

