

兔抗 MED1 多克隆抗体

中文名称：兔抗 MED1 多克隆抗体

英文名称： Anti-MED1 rabbit polyclonal antibody

别名： mediator complex subunit 1

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： MED1

宿主： Rabbit

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	The activation of gene transcription is a multistep process that is triggered by factors that recognize transcriptional enhancer sites in DNA. These factors work with co-activators to direct transcriptional initiation by the RNA polymerase II apparatus. The protein encoded by this gene is a subunit of the CRSP (cofactor required for SP1 activation) complex, which, along with TFIID, is required for efficient activation by SP1.
Applications:	WB, IHC
Name of antibody:	MED1
Immunogen:	Synthesized peptide derived from internal of human MED1.

Full name:	mediator complex subunit 1
Synonyms :	mediator complex subunit 1
SwissProt:	Q15648
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	168 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000



