

兔抗 HMG20B 多克隆抗体

中文名称：兔抗 HMG20B 多克隆抗体

英文名称：Anti-HMG20B rabbit polyclonal antibody

别名：SOXL; HMGX2; BRAF25; BRAF35; HMGXB2; PP7706; pp8857; SMARCE1r

抗原：HMG20B

储存：冷冻（-20℃）避光

宿主：Rabbit

反应种属：Human

相关类别：一抗

标记物：Unconjugate

克隆类型：Unconjugate

技术规格

Background:	Required for correct progression through G2 phase of the cell cycle and entry into mitosis. Required for RCOR1/C oREST mediated repression of neuronal specific gene promoters.
Applications:	WB
Name of antibody:	HMG20B
Immunogen:	Synthesized peptide derived from N-terminal of human HMG20B.
Full name:	high mobility group 20B
Synonyms :	SOXL; HMGX2; BRAF25; BRAF35; HMGXB2; PP7706; pp8857; SMARCE1r

SwissProt:	Q9P0W2
WB Predicted band size:	36 kDa
WB Positive control:	Jurkat cells lysate
WB Recommended dilution:	500-3000

