

## DTX3 抗原（重组蛋白）

中文名称：DTX3 抗原（重组蛋白）

英文名称：DTX3 Antigen (Recombinant Protein)

储存：冷冻（-20℃）

相关类别：抗原

概述

Fusion protein corresponding to C terminal 248 amino acids of human DTX3

技术规格

<b>Full name:</b>	deltex E3 ubiquitin ligase 3
<b>Synonyms:</b>	RNF154; deltex3
<b>Swissprot:</b>	Q8N9I9
<b>Gene Accession:</b>	BC114498
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	Regulator of Notch signaling, a signaling pathway involved in cell-cell communications that regulates a broad spectrum of cell-fate determinations. Probably acts both as a positive and negative regulator of Notch, depending on the developmental and cell context (By similarity). Functions as an ubiquitin ligase protein in vitro, suggesting that it may regulate the Notch pathway via some ubiquitin ligase activity.