

## 兔抗 MED23 多克隆抗体

中文名称：兔抗 MED23 多克隆抗体

英文名称：Anti-MED23 rabbit polyclonal antibody

别名：mediator complex subunit 23 SUR2; CRSP3; MRT18; SUR-2; ARC130; CRSP130; CRSP133; DRIP130

相关类别：一抗

储存：冷冻 (-20℃)

抗原：MED23

宿主：Rabbit

反应种属：Human, Mouse, Rat

标记物：Unconjugate


克隆类型：rabbit polyclonal

### 技术规格

**Background:**

The activation of gene transcription is a multistep process that is triggered by factors that recognize transcriptional enhancer sites in DNA. These factors work with co-activators to direct transcriptional initiation by the RNA polymerase II apparatus. The protein encoded by this gene is a subunit of the CRSP (cofactor required for SP1 activation) complex, which, along with TFIID, is required for efficient activation by SP1. This protein is also a component of other multisubunit complexes e.g. thyroid hormone receptor (TR-) associated proteins which interact with TR and facili

	itate TR function on DNA templates in conjunction with initiation factors and cofactors. This protein also acts as a metastasis suppressor. Several alternatively spliced transcript variants encoding different isoforms have been described for this gene.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	MED23
<b>Immunogen:</b>	Synthetic peptide of human MED23
<b>Full name:</b>	mediator complex subunit 23
<b>Synonyms:</b>	SUR2; CRSP3; MRT18; SUR-2; ARC130; CRSP130; CRSP133; DRIP130
<b>SwissProt:</b>	Q9ULK4
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human ovarian cancer;
<b>IHC Recommend dilution:</b>	30-150



**mlbio** 梅联生物  
Good elisakit producers