

兔抗 GAB2(Phospho-Ser623) 多克隆抗体

中文名称：兔抗 GAB2(Phospho-Ser623) 多克隆抗体

英文名称： Anti-GAB2(Phospho-Ser623) rabbit polyclonal antibody

别名： GAB2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： GAB2(Phospho-Ser623)

反应种属： Human Mouse

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

This gene is a member of the GRB2-associated binding protein (GAB) gene family. These proteins contain pleckstrin homology (PH) domain, and bind SHP2 tyrosine phosphatase and GRB2 adapter protein. They act as adapters for transmitting various signals in response to stimuli through cytokine and growth factor receptors, and T- and B-cell antigen receptors. The protein encoded by this gene is the principal activator of phosphatidylinositol-3 kinase in response to activation of the high affinity IgE receptor. Two alternatively spliced transcripts encoding different isoforms have b

	een described for this gene. [provided by RefSeq, No v 2009]
Applications:	WB, IHC
Name of antibody:	GAB2(Phospho-Ser623)
Immunogen:	Peptide sequence around phosphorylation site of Serine 623(P-S-S(p)-P-S) derived from Human Gab2 .
Full name:	GRB2-associated binding protein 2
Synonyms :	GAB2
SwissProt:	Q9UQC2
IHC positive control:	Human breast carcinoma tissue
IHC Recommend dilution:	50-100
WB Predicted band size:	74 kDa
WB Positive control:	293 cells, Hela cells and HepG2 cells lysates
WB Recommended dilution:	500-1000



