

兔抗 HMGCLL1 多克隆抗体

中文名称：兔抗 HMGCLL1 多克隆抗体

英文名称：Anti-HMGCLL1 rabbit polyclonal antibody

别名：3-hydroxymethyl-3-methylglutaryl-CoA lyase like 1; er-cHL; bA418P12.1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：HMGCLL1

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Non-mitochondrial 3-hydroxymethyl-3-methylglutaryl-CoA lyase that catalyzes a cation-dependent cleavage of (S)-3-hydroxy-3-methylglutaryl-CoA into acetyl-CoA and acetoacetate, a key step in ketogenesis, the products of which support energy production in nonhepatic animal tissues.
Applications:	ELISA, IHC
Name of antibody:	HMGCLL1
Immunogen:	Fusion protein of human HMGCLL1
Full name:	3-hydroxymethyl-3-methylglutaryl-CoA lyase like 1
Synonyms:	er-cHL; bA418P12.1
SwissProt:	Q8TB92

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human colorectal cancer
IHC Recommend dilution:	50-300



