

兔抗 MED6 多克隆抗体

中文名称：兔抗 MED6 多克隆抗体

英文名称： Anti-MED6 rabbit polyclonal antibody

别名： ARC33; NY-REN-28

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： MED6

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Mediator of RNA polymerase II transcription subunit 6 is one of the subunits of the Mediator complex. It is an enzyme that in humans is encoded by the MED6 gene. Component of the Mediator complex, a coactivator involved in the regulated transcription of nearly all RNA polymerase II-dependent genes. Mediator functions as a bridge to convey information from gene-specific regulatory proteins to the basal RNA polymerase II transcription machinery.
Applications:	ELISA, WB, IHC
Name of antibody:	MED6
Immunogen:	Full length fusion protein

Full name:	mediator complex subunit 6
Synonyms :	ARC33; NY-REN-28
SwissProt:	O75586
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and human lung cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	28 kDa
WB Positive control:	Hela cell
WB Recommended dilution:	200-1000



